MEDOVATIONS TRADITIONAL™ MERCURY-FILLED MALONEY & HURST BOUGIES

For technical information and/or
information regarding safety and
effectiveness, please call 800-558-6408

Caution: Federal (U.S.A.) Law restricts this device to
sale by or on the order of a physician.

102 E Keefe Ave
Milwaukee, WI 53212 USA
414.265.7620

Fax: 414.265.7628

e-mail: medo@medovations.com

1.800.558.6408
www.medovations.com

Made in the USA

MEDOVATIONS TRADITIONAL™ MERCURY-FILLED MALONEY & HURST BOUGIES

*Directions for Use and Cleaning.
Read carefully prior to use.*

DESCRIPTION

TRADITIONAL™ MERCURY-FILLED MALONEY & HURST BOUGIES

This package contains one or more polyvinylchloride (PVC) Bougies. Maloney Bougies have a tapered distal end. Hurst Bougies have a blunt distal end.

Traditional™ Bougies are available in 23 sizes from 16 French to 60 French, in 2 French increments (e.g. 16FR, 18FR, 20FR...). All sizes can be purchased individually, as a complete set (contains all 23 sizes), or as a mini –set (contains 10 sizes).

INTENDED USE

For dilatation of upper esophageal webs, lower esophageal rings, caustic strictures, peptic esophageal strictures, and temporary ease of esophageal carcinoma.

INSPECTION – PRIOR TO USE

Carefully inspect each Traditional™ Bougie prior to use:

1. Inspect new Bougies for damage from shipping and handling.
2. Inspect older Bougies for wear. To do this, bend and rotate the bougie at the distal tip, shaft and proximal end. Worn bougies should NOT be used. Signs of wear can include: surface cracks, cuts, and extremely faded print.
3. Make sure the three (3) year expiration date has not been exceeded.

CONTRAINDICATIONS

Contraindications include those specific to upper GI endoscopy.

Contraindications to dilation include, but are not limited to: uncooperative patient; asymptomatic strictures; inability to advance the bougie through the strictured area; coagulopathy; known or suspected perforation; severe inflammation or scarring near the dilation site, recent myocardial infarction, active ulcer and severe cervical arthritis.

POTENTIAL COMPLICATIONS

Potential complications associated with upper gastrointestinal endoscopy and esophageal dilation include, but are not limited to: perforation, hemorrhage, aspiration, fever, infection, allergic reaction to medication, hypotension, respiratory depression or arrest, cardiac arrhythmia or arrest.

INSTRUCTIONS FOR USE

It is recommended that this procedure be performed under fluoroscopy. Traditional™ Bougies are completely radiopaque.

1. Perform screening endoscopy to identify strictured area.
2. Generously lubricate bougie and advance it to strictured area.
 - a. **Note:** Medovations Traditional™ Bougies have two (2) series* of numbers printed on each bougie on opposing sides. One series measures the distance to the distal tip, identified as **FROM TIP (American)**. The other series measures the distance to the stricture, identified as **FROM MAX O.D. (European)**. These printed distance markings can be used as a guide to avoid over-extending into the gastric cavity.
3. Proceed with esophageal dilation. Determine the bougie to be used by the size of the strictures, starting with a smaller French sized bougie and gradually increasing the French size. Hold the bougie at the blunt, proximal end.
4. Upon completion of esophageal dilation, remove bougie from patient.

Note: Successive dilations can be accomplished by repeating steps 2-4.

5. Clean and prepare bougie for future use following cleaning instructions below.

*Note: The Hurst style bougie has only one series of printed numbers, FROM TIP (American).

CLEANING INSTRUCTIONS

CAUTION:

1. During cleaning, inspect integrity and function of device to determine advisability of reuse. If kinks, excessive wear, bends or breaks exist, do not use. Dispose of device per institutional guidelines for biohazardous medical waste.
2. **DO NOT USE** Betadine to clean the bougies as it may cause discoloration of the material.
3. Traditional™ Bougies **CANNOT** be steam autoclaved. The high temperatures of steam autoclaving will destroy the bougie.
4. **DO NOT USE** alcohol for cleaning or disinfection.
5. Appropriate gloves, gowns, and other protective wear should be worn when cleaning and sterilizing soiled or used Bougies.
6. It is the responsibility of the user to qualify any deviations from the recommended methods of cleaning.

Bougies must be thoroughly cleaned prior to use.

PRE-CLEANING:

1. Immediately after use, soak bougie in a solution of warm water and enzymatic cleaning solution. **Note:** Refer to the enzymatic cleaning solution instructions for recommended duration of cleaning.
2. Remove gross contaminants on the external surfaces of the bougie by hand, using a soft bristle brush, sponge, or wipe, with warm water and enzymatic solution.
3. Thoroughly rinse the bougie with clean water.
4. Pat dry with clean, non-contaminating wipes.
5. Proceed with high-level disinfection / sterilization of bougie.

HIGH-LEVEL DISINFECTION / STERILIZATION:

The following disinfection solutions can be used with the Medovations Traditional™ Bougies:

- Glutaraldehyde
- Ortho-phthalaldehyde (OPA)
- Peracetic acid
- Hydrogen peroxide

Please refer to the manufacturer's instructions for complete information and precautions for disinfection solutions with any of the above active ingredients.

The Bougies have additionally been validated with the following specific brands:

- Cidex® OPA High Level Disinfectant
- Steris Resert™ XL HLD High Level Disinfectant
- Rapicide PA High Level Disinfectant
- Reliance™ EPS Endoscope Processing System
- Medivators Advantage Plus System
- Medivators DSD series of automated endoscope reprocessors
- Steris SYSTEM 1E® Liquid Chemical Sterilant Processing System
- STERRAD NX® System. The STERRAD NX® sterilization validation is based on a standard 28 minute cycle. The maximum load limit for the STERRAD NX® System is 4850 grams (10.7 lbs) per shelf. Do not use sterilization pouches/peel pouches when performing the STERRAD NX®

sterilization process. Doing so may result in damage to the device. It is recommended to use polypropylene sterile wrapped trays.

Please refer to the manufacturer's instructions for complete information and precautions.

Note: Proper rinsing with clean water and complete drying of the bougie is necessary with any disinfection solution.

STORAGE

The Traditional™ Bougies should be stored flat, out of direct light, and at room temperature. The Medovations Bougie Storage Case provides a protected environment.

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