

CareGuard™

A bold new look
for the common bite block



medovations, inc.
102 east keefe avenue, milwaukee, wi 53212 usa
800.558.6408, 414.265.7620, fax 414.265.7628
medo@medovations.com
www.medovations.com

Why compromise sensitive, expensive equipment when a simple solution from Medovations can change everything? Our new line of bite blocks takes a colorful twist on caring for your endoscopes and your patients.

At a glance you'll know it's different

Medovations' CareGuard™ Bite Blocks come in two fresh, color-coded sizes—48 FR and 60 FR. You'll reach for the right size every time.

Durable, dependable design

Created with the patient and the caregiver in mind, CareGuard Bite Blocks offer a variety of features. An integrated dental retention rim maintains the position of the bite block during procedures, and lateral ports on the 60 FR bite block allow for additional suctioning. Choose from optional, adjustable straps that securely and comfortably fit any patient.

The 60 FR CareGuard Bite Block easily accommodates a 60 FR Medovations Esophageal Bougie or Wire Guided Dilator. Trust the construction of CareGuard Bite Blocks to provide the protection you need.

Disposable and latex-free

Individually wrapped CareGuard Bite Blocks are made for one-time use, reducing the risk of cross contamination.



You know Medovations for our innovative, reliable, and cost effective solutions. Now, there's another reason to smile.

Part No.	Size	Description	Qty/Case
CareGuard™ Bite Block			
1197-48	48 FR	With dental rim	100
1197-60	60 FR	With dental rim and lateral ports	100
CareGuard™ Bite Block with Latex-free Blue Strap			
1196-48	48 FR	With dental rim	100
1196-60	60 FR	With dental rim and lateral ports	100
CareGuard™ Bite Block with Latex-free Soft Strap			
1198-48	48 FR	With dental rim	100
1198-60	60 FR	With dental rim and lateral ports	100

Made in the USA



48 FR




60 FR

CareGuard™

The bold new choice for upper GI procedures

To learn more, visit medovations.com or call 800-558-6408 to speak with a personal representative.



Dental retention rim

Lateral ports offer additional access