

Solutions at a price that
won't clean out your budget.



You know Medovations for its broad lines of gastroenterology and surgical products. Our new line of CleanFreak™ brushes is designed to meet the cleaning requirements of flexible endoscopes:

- Range of brush sizes and configurations
- 1.7 mm catheter for control and flexibility
- Single-use durable design
- Non-sterile and Latex-free
- Made in the USA



If it has to be clean,
leave it to a CleanFreak.

Clean
freak™



Clean
freak™

 medovations®

medovations, inc.
102 east keefe avenue, milwaukee, wi 53212 usa
800.558.6408, 414.265.7620, fax 414.265.7628
medo@medovations.com
www.medovations.com

 medovations®

End to end,
bristle for bristle,
there's nothing
quite like
CleanFreak.™

What do you expect from your endoscope cleaning brushes? How about a broad line from one supplier, and confidence in the products that help you clean sensitive, expensive equipment. CleanFreak Endoscope Cleaning Brushes are precisely what you expect from Medovations—solutions at a price that won't clean out your budget. End to end, bristle for bristle, you'll feel better with a CleanFreak by your side.

Medovations' tradition of innovative endoscopic devices continues with CleanFreak Endoscope Cleaning Brushes. To learn more, contact a personal product representative at 800.558.6408.

Standard Brush | 1200-01

This single-head brush is a must have accessory, combining a non-slip catheter with a tapered head and a protective tip. If one brush ever did the work of many, this would be it.

Extra Wide Brush | 1200-02

This brush contains longer bristles, a plastic tip, and a compact head that provides optimal cleaning capability for channel diameters between 2.8 mm to 6 mm.

Extra Long Brush | 1200-03

The tapered brush head and firm yellow catheter flexes easily when guided through the channel. With a longer head design, this brush doubles the cleaning contact with every pass through the channel.

Dual Ended Brush | 1200-04

This dual-ended brush enables efficient cleaning by using two brushes in one pass. The textured catheter helps maintain a secure grip.

Combination Brush | 1200-05

For maximum convenience, a standard brush and a valve brush are now in one product. A bright yellow tip on the channel brush helps protect the channel from abrasion.

Valve Brush | 1200-06

This compact cleaning brush features a densely packed head of soft, gentle bristles to clean the delicate areas around the control head and valve cavities. The handle design ensures a comfortable, secure grip.

CleanFreak™ Endoscope Brushes for Flexible Endoscopes

Part No.	Description	Bristle Diameter	Sheath Diameter	Length	Minimum Channel Dia	Qty/Case
1200-01	Standard Brush	5–7 mm	1.7 mm	230 cm	2.0 mm	50
1200-02	Extra Wide Brush	7–9 mm	1.7 mm	230 cm	2.8 mm	50
1200-03	Extra Long Brush	5–7 mm	1.7 mm	230 cm	2.0 mm	50

Single-use, non-sterile. Made in the USA.

CleanFreak™ Endoscope Brushes for Flexible Endoscopes

Part No.	Description	Bristle Diameter	Sheath Diameter	Length	Minimum Channel Dia	Qty/Case
1200-04	Dual Ended Brush	5–7 mm	1.7 mm	230 cm	2.0 mm	50
1200-05	Combination Brush (Valve & Channel Brush)	5–7 mm	1.7 mm	230 cm	2.0 mm	50
1200-06	Valve Brush	1.0 cm				50

Single-use, non-sterile. Made in the USA.