



medovations

## Esophageal Dilatation Systems



medovations

medovations, inc.  
102 east keefe avenue, milwaukee, wi 53212 usa  
800.558.6408, 414.265.7620, fax 414.265.7628  
medo@medovations.com  
[www.medovations.com](http://www.medovations.com)



Medovations Innovative products from a new perspective


It's simple. At Medovations, we have a way of bringing the right product to the marketplace because we think about your needs. All of which is evidenced by the little things we do. The results? An unsurpassed range of products and services. Whether it's a new device or a subtle enhancement, there's always a better way when innovation is part of your name.

**M-Flex™ Blue Silicone Bougies**  
(Tungsten-filled)

FR Size (mm)	Maloney (tapered) Tip Part No.	Hurst (blunt) Tip Part No.
20	1206-20	1207-20
22	1206-22	1207-22
24	1206-24	1207-24
26	1206-26	1207-26
28	1206-28	1207-28
30	1206-30	1207-30
32	1206-32	1207-32
34	1206-34	1207-34
36	1206-36	1207-36
38	1206-38	1207-38
40	1206-40	1207-40
42	1206-42	1207-42
44	1206-44	1207-44
46	1206-46	1207-46
48	1206-48	1207-48
50	1206-50	1207-50
52	1206-52	1207-52
54	1206-54	1207-54
56	1206-56	1207-56
58	1206-58	1207-58
60	1206-60	1207-60

**M-Flex™ Dilatation Systems**

Set	Part No.	Contents
Maloney (tapered) 21-piece full set	1206-00	21 bougies, 20–60 FR
Maloney (tapered) 10-piece mini-set	1206-10	10 bougies, 36–54 FR
Hurst (blunt) 21-piece full set	1207-00	21 bougies, 20–60 FR
Hurst (blunt) 10-piece mini-set	1207-10	10 bougies, 36–54 FR
Bougie Storage Case	1210-01	Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H



**M-Flex™ Blue Silicone Bougies**



Store and transport M-Flex™ bougies in our storage case, which holds up to two full sets. A plexiglass window lets you identify the FR size on the end of each bougie, or remove the window to access every product without opening the case.



**M-Flex™ Ultra-flexible Silicone Bougies.**

Our silicone bougies feature added flexibility that some users prefer. Tungsten-filled and latex-free, there's nothing quite like the M-Flex™.

**Benefits**

- Ultra-flexible tapered tip (Maloney)
- Completely radiopaque from distal tip to proximal end
- American and European depth markings
- Mercury-free and Latex-free

**Features**

- 21 sizes: 20 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini-set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No-charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)

# SafeGuide™ Wire Guided Dilators

## SafeGuide™ Over the Guidewire Esophageal Dilators

FR Size (mm)	Part No.	
15	5	1214-15
18	6	1214-18
21	7	1214-21
24	8	1214-24
27	9	1214-27
30	10	1214-30
33	11	1214-33
36	12	1214-36
39	13	1214-39
42	14	1214-42
45	15	1214-45
48	16	1214-48
51	17	1214-51
54	18	1214-54
57	19	1214-57
60	20	1214-60

## SafeGuide™ Dilatation Systems

Set	Part No.	Contents
16-piece full set	1214-00	16 dilators, 15-60 FR, 1 storage case, 1 marked spring tip guidewire, 1 cleaning brush, 25 disposable cleaning adapters
7-piece mini-set	1214-10	7 dilators (15, 21, 27, 33, 39, 42, 45 FR), 1 storage case, 1 marked spring tip guidewire, 1 cleaning brush, 25 disposable cleaning adapters

## SafeGuide™ Accessories

Part	Part No.	Contents
SafeGuide™ Storage Case	1214-01	Heavy-duty storage case for up to 16 SafeGuide™ Dilators. 36"L x 10"W x 4"H
Marked Spring Tip Guidewires	1214-02	Box of 2 guidewires for use with all dilators. 210cm
Cleaning Brushes	1214-03	Box of 3 cleaning brushes to pass through inner lumen of all dilators
Disposable Cleaning Adapters	1214-04	Box of 25 disposable cleaning adapters to flush inner lumen of all dilators

### SafeGuide™ Flexibility and control for exacting procedures.

At times, there is need for greater accuracy. That's why there's SafeGuide™, a remarkably flexible over the guidewire solution that ensures exact placement of the dilator. Position the guidewire with the assistance of an endoscope and direct the SafeGuide™ dilator where you need it.

#### Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Latex-free
- SafeGuide™ guidewire features radiopaque bands

#### Features

- 16 sizes: 15 FR to 60 FR (in 3 FR increments)
- Choose a full set, mini-set or individual pieces
- Easy to clean with optional brushes and adapters
- With proper care, reusable SafeGuide™ dilators provide years of service

The single-piece construction of the SafeGuide™ Spring Tip Guidewire resists breakage and flexes easily

Optional Marked Spring Tip Guidewire, Cleaning Brush, Cleaning Adapters and Storage Case are available



Store and transport SafeGuide™ Dilators in a convenient storage case. A labeled slot holds each dilator secure. The case holds one full set, 16 dilators total.

**WeightRight™ Mercury-free Bougies**  
(Tungsten-filled)

FR Size (mm)	Maloney (tapered) Tip Part No.	Hurst (blunt) Tip Part No.
16	1208-16	1209-16
18	1208-18	1209-18
20	1208-20	1209-20
22	1208-22	1209-22
24	1208-24	1209-24
26	1208-26	1209-26
28	1208-28	1209-28
30	1208-30	1209-30
32	1208-32	1209-32
34	1208-34	1209-34
36	1208-36	1209-36
38	1208-38	1209-38
40	1208-40	1209-40
42	1208-42	1209-42
44	1208-44	1209-44
46	1208-46	1209-46
48	1208-48	1209-48
50	1208-50	1209-50
52	1208-52	1209-52
54	1208-54	1209-54
56	1208-56	1209-56
58	1208-58	1209-58
60	1208-60	1209-60

**WeightRight™ Dilatation Systems**

Set	Part No.	Contents
Maloney (tapered) 23-piece full set	1208-00	23 bougies, 16–60 FR
Maloney (tapered) 10-piece mini-set	1208-10	10 bougies, 36–54 FR
Hurst (blunt) 23-piece full set	1209-00	23 bougies, 16–60 FR
Hurst (blunt) 10-piece mini-set	1209-10	10 bougies, 36–54 FR
Bougie Storage Case	1210-01	Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H

# WeightRight™ Mercury-free Bougies



Store and transport WeightRight™ bougies in our storage case. It holds up to two full sets, 46 bougies total. A plexiglass window lets you identify the FR size on the end of each bougie, or remove the window to access every product without opening the case.



**WeightRight™. Heavy-weight control and convenience combined.**

Our tungsten-filled bougie mimics the weight and feel of mercury-filled bougies without the concerns of mercury usage and disposal. It's the weight you expect, and the convenience you demand.


**Benefits**

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Latex-free

**Features**

- 23 sizes: 16 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini-set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No-charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)

# Traditional™ Mercury-filled Bougies



### Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Latex-free

### Features

- 23 sizes: 16 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini-set or individual pieces
- Three-year warranty on product and workmanship
- No-charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)

### Traditional™. For some, there's no substitute.

Our classic mercury-filled bougie is still a favorite among many who prefer a mercury-weighted product. Our trade-in program makes disposal safe and easy. Yes, it's traditional by name, but there's nothing traditional about it.

### Traditional™ Mercury-filled Bougies

FR Size (mm)	Maloney (tapered) Tip Part No.	Hurst (blunt) Tip Part No.	
16	5.3	1210-16	1212-16
18	6	1210-18	1212-18
20	6.6	1210-20	1212-20
22	7.4	1210-22	1212-22
24	8	1210-24	1212-24
26	8.6	1210-26	1212-26
28	9.3	1210-28	1212-28
30	10	1210-30	1212-30
32	10.7	1210-32	1212-32
34	11.2	1210-34	1212-34
36	12	1210-36	1212-36
38	12.5	1210-38	1212-38
40	13	1210-40	1212-40
42	14	1210-42	1212-42
44	14.5	1210-44	1212-44
46	15.2	1210-46	1212-46
48	16	1210-48	1212-48
50	16.5	1210-50	1212-50
52	17.2	1210-52	1212-52
54	18	1210-54	1212-54
56	18.5	1210-56	1212-56
58	19.2	1210-58	1212-58
60	20	1210-60	1212-60

### Traditional™ Dilatation Systems

Set	Part No.	Contents
Maloney (tapered) 23-piece full set	1210-00	23 bougies, 16–60 FR
Maloney (tapered) 10-piece mini-set	1210-10	10 bougies, 36–54 FR
Hurst (blunt) 23-piece full set	1212-00	23 bougies, 16–60 FR
Hurst (blunt) 10-piece mini-set	1212-10	10 bougies, 36–54 FR
Bougie Storage Case	1210-01	Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H



With room for 46 bougies, the storage case for Traditional™ bougies ensures proper storage and easy transport. The FR size can be identified on the end of the bougie through a plexiglass window, or remove the window to access every bougie without opening the case.