



medovations®

# Esophageal Dilatation Systems



M-Flex™  
Ultra-flexible Silicone Bougies.



Store and transport M-Flex™ bougies in our storage case which holds up to two full sets. A plexiglass window enables you to identify the FR size on the end of each bougie, or remove the window to access every product without opening the case.

# M-Flex™ Blue Silicone Bougies.

Our silicone bougies feature added flexibility that some users prefer. Tungsten-filled and latex-free, there's nothing quite like the M-Flex™.

## Benefits

- Ultra-flexible tapered tip (Maloney)
- Completely radiopaque from distal tip to proximal end
- American and European depth markings
- Mercury-free and Latex-free

## Features

- 21 sizes: 20 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini-set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No-charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)

## M-Flex™ Dilatation Systems

Set	Part No.	Contents
Maloney (tapered) 21-piece full set	1206-00	21 bougies, 20–60 FR
Maloney (tapered) 10-piece mini-set	1206-10	10 bougies, 36–54 FR
Hurst (blunt) 21-piece full set	1207-00	21 bougies, 20–60 FR
Hurst (blunt) 10-piece mini-set	1207-10	10 bougies, 36–54 FR
Bougie Storage Case	1210-01	Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H

## M-Flex™ Blue Silicone Bougies

(Tungsten-filled)		Maloney (tapered)	Hurst (blunt)
FR Size	(mm)	Tip Part No.	Tip Part No.
20	6.6	1206-20	1207-20
22	7.4	1206-22	1207-22
24	8	1206-24	1207-24
26	8.6	1206-26	1207-26
28	9.3	1206-28	1207-28
30	10	1206-30	1207-30
32	10.7	1206-32	1207-32
34	11.2	1206-34	1207-34
36	12	1206-36	1207-36
38	12.5	1206-38	1207-38
40	13	1206-40	1207-40
42	14	1206-42	1207-42
44	14.5	1206-44	1207-44
46	15.2	1206-46	1207-46
48	16	1206-48	1207-48
50	16.5	1206-50	1207-50
52	17.2	1206-52	1207-52
54	18	1206-54	1207-54
56	18.5	1206-56	1207-56
58	19.2	1206-58	1207-58
60	20	1206-60	1207-60



# WeightRight™ Mercury-free Bougies.

Our tungsten-filled bougie mimics the weight and feel of mercury-filled bougies without the concerns of mercury usage and disposal. It's the weight you expect, and the convenience you demand.

## Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Latex-free

## Features

- 23 sizes: 16 FR to 60 FR (in 2 FR increments)
- Choose a full set, mini-set or individual pieces
- Three-year warranty on product and workmanship
- Easy disposal
- No-charge trade-in program for proper disposal of expired mercury-filled bougies (with purchase of new product)

## WeightRight™ Dilatation Systems

Set	Part No.	Contents
Maloney (tapered) 23-piece full set	1208-00	23 bougies, 16–60 FR
Maloney (tapered) 10-piece mini-set	1208-10	10 bougies, 36–54 FR
Hurst (blunt) 23-piece full set	1209-00	23 bougies, 16–60 FR
Hurst (blunt) 10-piece mini-set	1209-10	10 bougies, 36–54 FR
Bougie Storage Case	1210-01	Heavy-duty storage for up to 46 bougies. 32"L x 7"W x 4"H

## WeightRight™ Mercury-free Bougies

(Tungsten-filled)

FR Size	(mm)	Maloney (tapered) Tip Part No.	Hurst (blunt) Tip Part No.
16	5.3	1208-16	1209-16
18	6	1208-18	1209-18
20	6.6	1208-20	1209-20
22	7.4	1208-22	1209-22
24	8	1208-24	1209-24
26	8.6	1208-26	1209-26
28	9.3	1208-28	1209-28
30	10	1208-30	1209-30
32	10.7	1208-32	1209-32
34	11.2	1208-34	1209-34
36	12	1208-36	1209-36
38	12.5	1208-38	1209-38
40	13	1208-40	1209-40
42	14	1208-42	1209-42
44	14.5	1208-44	1209-44
46	15.2	1208-46	1209-46
48	16	1208-48	1209-48
50	16.5	1208-50	1209-50
52	17.2	1208-52	1209-52
54	18	1208-54	1209-54
56	18.5	1208-56	1209-56
58	19.2	1208-58	1209-58
60	20	1208-60	1209-60



Traditional™ Mercury-filled Bougies are still available in sets and as individual pieces. Contact us at 800.558.6408 to order.



WeightRight™  
Heavy-weight control and convenience combined.



Store and transport WeightRight™ bougies in our storage case. It holds up to two full sets, 46 bougies total. A plexiglass window enables you to identify the FR size on the end of each bougie, or remove the window to access every product without opening the case.



SafeGuide™  
Flexibility and control for exacting procedures.



Store and transport SafeGuide™ dilators in a convenient storage case. A labeled slot holds each dilator secure. The case holds one set, 16 dilators total.

# SafeGuide™ Wire Guided Dilators.

At times, there is need for greater accuracy. That's why there's SafeGuide™ a remarkably flexible over the guidewire solution that ensures exact placement of the dilator. Position the guidewire with the assistance of an endoscope and direct the SafeGuide™ dilator where you need it.

## Benefits

- Completely radiopaque from distal tip to proximal end
- Flexible material does not attract particulate
- American and European depth markings
- Latex-free
- SafeGuide™ guidewire features radiopaque bands

## Features

- 16 sizes: 15 FR to 60 FR (in 3 FR increments)
- Choose a full set, mini-set or individual pieces
- Easy to clean with optional brushes and adapters
- With proper care, reusable SafeGuide™ dilators provide years of service

## SafeGuide™ Dilatation Systems


Set	Part No.	Contents
16-piece full set	1214-00	16 dilators, 15–60 FR, 1 storage case, 1 marked spring tip guidewire, 1 cleaning brush, 25 disposable cleaning adapters
7-piece mini-set	1214-10	7 dilators (15, 21, 27, 33, 39, 42, 45 FR), 1 storage case, 1 marked spring tip guidewire, 1 cleaning brush, 25 disposable cleaning adapters

## SafeGuide™ Accessories

Part	Part No.	Contents
SafeGuide™ Storage Case	1214-01	Heavy-duty storage case for up to 16 SafeGuide™ Dilators. 36"L x 10"W x 4"H
Marked Spring Tip Guidewires	1214-02	Box of 2 guidewires for use with all dilators. 210cm
Reusable Cleaning Brushes	1214-03	Box of 3 cleaning brushes to pass through inner lumen of all dilators
Disposable Cleaning Adapters	1214-04	Box of 25 disposable cleaning adapters to flush inner lumen of all dilators

## SafeGuide™ Over the Guidewire Esophageal Dilators

FR Size	(mm)	Part No.
15	5	1214-15
18	6	1214-18
21	7	1214-21
24	8	1214-24
27	9	1214-27
30	10	1214-30
33	11	1214-33
36	12	1214-36
39	13	1214-39
42	14	1214-42
45	15	1214-45
48	16	1214-48
51	17	1214-51
54	18	1214-54
57	19	1214-57
60	20	1214-60



The single-piece construction of the SafeGuide™ Spring Tip Guidewire resists breakage and flexes easily.

Optional Marked Spring Tip Guidewire, Cleaning Brush, Cleaning Adapters and Storage Case are available.



medovations®

medovations, inc.  
102 east keefe avenue, milwaukee, wi 53212 usa  
800.558.6408, 414.265.7620, fax 414.265.7628  
medo@medovations.com  
[www.medovations.com](http://www.medovations.com)

M-Flex™ SafeGuide™ WeightRight™ and Traditional™ are trademarks of Medovations, Inc.  
©2010 Medovations, Inc. All rights reserved. 080-0021 7/10