

Aspiration Tubing



Aspiration Tubing

High-flow, latex-free aspiration tubing withstands the powerful demands of Power Assisted Lipoplasty. Ten feet of proprietary, "memory-free" tubing provides for maximum freedom of movement.



medovations

Model	Sterile Tubing Description	Count Per Case
SU1159	Aspiration Tubing, Tapered MicroAire® Compatible, sterile	10
SU1160	Aspiration Tubing, sterile	10



SU1159



SU1160

Endoscopy • Liposuction • Neurosurgery



Infiltration Tubing



Infiltration Tubing

Model SU1161—Dual Spike Tubing

Model SU1163—Single Spike Tubing

Large bore “memory-free” tubing in full eight foot lengths with a luer lock fitting on the operating end. Available in dual and single spike configurations.

Infiltration Tubing for Peristaltic Pumps

Model SU1168

Model SU1170-LY—Lysonix® compatible

Infiltration infusion pump tubing has a nine-foot length and a luer lock fitting. Peristaltic silicone tubing is pressure-locked to a six-foot section of tubing with a bag spike. The silicone tubing section is designed for increased flow and easy insertion on the pump head.



medovations®

Model	Description	Count Per Case
SU1161	Infiltration Dual Spike Tubing, sterile	10
SU1163	Infiltration Single Spike Tubing, sterile	10
SU1168	Infiltration Tubing for Peristaltic Pump, sterile	10
SU1170-LY	Bifurcated Infiltration Tubing for LySonix® Irrigator, sterile	10



SU1161



SU1163



SU1168



SU1170-LY

Endoscopy • Liposuction • Neurosurgery

