



Thoracic Catheter

- Straight
- Right Angle



- Clear, kink-resistant, medical grade PVC material
- Angle-cut end for easy insertion
- QuickConnect™ funnel end can be cropped to form integral funnel connector
- QuickConnect™ funnel connector provides quick, easy connection to standard drainage units
- Smooth drainage eyelets and open distal tip enhance efficient drainage and maximize patient comfort
- Radiopaque stripe
- Depth markings in 2 cm increments
- Right Angle style available
- Adult, pediatric, and infant sizes
- Sterile

Part Number	French Size	O.D. (mm)	Drainage Eyelets	Length	Qty per cs	Part Number	French Size	O.D. (mm)	Drainage Eyelets	Length	Qty per cs
STRAIGHT:						RIGHT ANGLE:					
1240-12	12	4.0	4	20"(51cm)	10	1245-16	16	5.3	4	20"(51cm)	10
1240-16	16	5.3	4	20"(51cm)	10	1245-20	20	6.7	5	20"(51cm)	10
1240-20	20	6.7	5	20"(51cm)	10	1245-24	24	8.0	6	20"(51cm)	10
1240-24	24	8.0	6	20"(51cm)	10	1245-28	28	9.3	6	20"(51cm)	10
1240-28	28	9.3	6	20"(51cm)	10	1245-32	32	10.7	6	20"(51cm)	10
1240-32	32	10.7	6	20"(51cm)	10	1245-36	36	12.0	6	20"(51cm)	10
1240-36	36	12.0	6	20"(51cm)	10	1245-40	40	13.3	6	20"(51cm)	10
1240-40	40	13.3	6	20"(51cm)	10						

102 E. KEEFE AVENUE • MILWAUKEE, WISCONSIN 53212 U.S.A.

Tel 414/265•7620 • Fax 414/265•7628

e-mail: medo@medovations.com

1•800•558•6408



Trocar Catheter



- Clear, kink-resistant, medical grade PVC material
- Polished, surgically sharp metal trocar with white hand piece
- Safety stop in distal end of cannula maintains proper position of metal trocar within the cannula and prevents trocar advancement
- QuickConnect™ funnel connector provides quick and easy connection to standard drainage units
- Radiopaque stripe
- Depth markings in 2cm increments
- Adult, pediatric, and infant sizes
- Sterile

Part Number	French Size	O.D. (mm)	Length	Qty per cs
1250-08	8	2.7	9"(23cm)	10
1250-10	10	3.3	9"(23cm)	10
1250-12	12	4.0	9"(23cm)	10
1250-16	16	5.3	10"(25cm)	10
1250-20	20	6.7	16"(41cm)	10
1250-24	24	8.0	16"(41cm)	10
1250-28	28	9.3	16"(41cm)	10
1250-32	32	10.7	16"(41cm)	10

102 E. KEEFE AVENUE • MILWAUKEE, WISCONSIN 53212 U.S.A.


Tel 414/265•7620 • Fax 414/265•7628

e-mail: medo@medovations.com

1•800•558•6408



Gastric Sump Tube



Blue side arm caps to 5-in-1 connector when not in use, reducing reflux accidents

- Double lumen PVC tube - a large lumen for suction and a small lumen for irrigation and venting
- Blue side arm extends into the smaller irrigation/vent lumen
- Blue side arm allows simultaneous irrigation of stomach without discontinuing suction and doubles as a vent lumen when irrigation is not taking place
- Unique double lumen design and vent feature enhances patient safety by eliminating risk of injury to stomach lining
- Funnel end with removable 5-in-1 connector easily connects to suction units
- Blue side arm caps to 5-in-1 connector when not in use, reducing reflux accidents
- Radiopaque stripe facilitates proper placement under x-ray
- Adult and pediatric sizes
- Sterile

Part Number	French Size	O.D. (mm)	Length	Qty per cs
1217-10	10	3.3	36"(91cm)	50
1217-12	12	4.0	48"(122cm)	50
1217-14	14	4.7	48"(122cm)	50
1217-16	16	5.3	48"(122cm)	50
1217-18	18	6.0	48"(122cm)	50

102 E. KEEFE AVENUE • MILWAUKEE, WISCONSIN 53212 U.S.A.


Tel 414/265•7620 • Fax 414/265•7628

e-mail: medo@medovations.com

1•800•558•6408



Flexible Yankauer Suction Device



- Flexible version of the traditional Yankauer suction device
- Suction tip can be formed to a precise suctioning angle during procedure
- Flexible suction end minimizes tissue trauma during suctioning
- Thermoformed suction end returns to original angle when not in use
- Available in 3 volume capacities; fine, regular, and high
- Suction Control Port (available on regular capacity) gives end-user complete control of vacuum pressure
- Basket Tip Guard (available on regular capacity) prevents tissue and bone chips from occluding the device
- Sterile

Suction Control Port and Basket Tip Guard options

Part Number	Suction Capacity	Suction Control Port	Basket Tip Guard	Qty per cs
1283-01	Fine	No	No	50
1284-01	Regular	No	No	50
1284-02	Regular	No	Yes	50
1284-03	Regular	Yes	No	50
1284-04	Regular	Yes	Yes	50
1285-01	High	No	No	50

102 E. KEEFE AVENUE • MILWAUKEE, WISCONSIN 53212 U.S.A.

Tel 414/265•7620 • Fax 414/265•7628

e-mail: medo@medovations.com

1•800•558•6408



Universal Funnel Tubing



- Clear, kink-resistant PVC
- "Bubble" every 3ft (1m) provides QuickConnect™ integral connector for devices, suction units, or tubing
- Packaged in 100ft (30m) roll in an easy-dispensing storage box
- High resistance to vacuum collapse



"Bubble" every 3ft (1m)



Integral connector formed when "Bubble" is cut


Part Number	I.D.	Length	Qty per cs
1268-01	3/16"(5mm)	100'(30m)	1 Roll
1268-02	1/4"(6mm)	100'(30m)	1 Roll
1268-03	9/32"(7mm)	100'(30m)	1 Roll

5-in-1 Connectors

- Stepped
- Smooth



STEPPED



- Inner diameter expands proportionately with outer diameter
- Polypropylene material can be easily cut to match ID of connector to ID of tube
- Cut connector maximizes flow rate and prevents fluid buildup at connection point
- Reduces inventory of connectors by accommodating five different tube sizes
- 5 connection points on 1 connector
- Sterile or non-sterile

SMOOTH



- Reduces inventory of connectors by accommodating five different tube sizes
- 5 connection points on 1 connector
- Non-sterile

Part Number	Description	Material	Fits Tubes	Qty per cs
<i>Stepped:</i>				
1350-00	Non-Sterile	Polypropylene	5/16" - 9/16" (8mm-14mm)	100
1350-55	Sterile	Polypropylene	5/16" - 9/16" (8mm-14mm)	100
<i>Smooth:</i>				
1350-75	Non-Sterile	Polypropylene	3/16" - 7/16" (5mm-11mm)	100

102 E. KEEFE AVENUE • MILWAUKEE, WISCONSIN 53212 U.S.A.

Tel 414/265•7620 • Fax 414/265•7628

e-mail: medo@medovations.com

1•800•558•6408