

Product Specifications

Cables for use with Circon-ACMI/ Olympus and other light sources

40 FR

Part No. 1450-40, InnerVision™ Fiber Optic Cable, 40 FR, Circon-ACMI/Olympus Connector. Clean, non-sterile.

50 FR

Part No. 1450-50, InnerVision™ Fiber Optic Cable, 50 FR, Circon-ACMI/Olympus Connector. Clean, non-sterile.

56 FR

Part No. 1450-56, InnerVision™ Fiber Optic Cable, 56 FR, Circon-ACMI/Olympus Connector. Clean, non-sterile.

Wolf/Dyonics Adapter for Circon-ACMI/Olympus light source

Part No. 1450-02, Wolf/Dyonics Adapter for InnerVision™ Cables with a Circon-ACMI/Olympus Connector. Clean, non-sterile.

Cables for use with Storz light source

40 FR

Part No. 1455-40, InnerVision™ Fiber Optic Cable, 40 FR, Storz Connector. Clean, non-sterile.

50 FR

Part No. 1455-50, InnerVision™ Fiber Optic Cable, 50 FR, Storz Connector. Clean, non-sterile.

56 FR

Part No. 1455-56, InnerVision™ Fiber Optic Cable, 56 FR, Storz Connector. Clean, non-sterile.

InnerVision™ Detachable Tips

40 FR

Part No. 1465-40, InnerVision™ Detachable Tips, 40 FR, box of 5, sterile.

50 FR

Part No. 1465-50, InnerVision™ Detachable Tips, 50 FR, box of 5, sterile.

56 FR

Part No. 1465-56, InnerVision™ Detachable Tips, 56 FR, box of 5, sterile.

Light Source Compatibility Guide

InnerVision™ Fiber Optic Cable with Circon-ACMI/Olympus Connector is compatible with:

- ACMI
- AFI
- ASPEN
- CIRCON
- CUDA
- EDER
- KLI
- LUXTEC
- OLYMPUS
- PILLING
- STRYKER
- V. MUELLER
- WEHMER
- ZIMMER

With a Wolf/Dyonics Adapter the InnerVision™ Fiber Optic Cable with Circon-ACMI/Olympus Connector is compatible with:

- WOLF
- DYONICS

InnerVision™ Fiber Optic Cable with Storz Connector is compatible with:

- KARL STORZ (5.0 Maximum Bundle Diameter)